

Boots & Cables | Programming connection



Product*	RIC R**	RIC 13	RIC 312	MICRO RIC 312	BTE 13 (Standard)	BTE 13 (Power Plus)	MINI BTE 312	Custom
Livio™ AI/ Livio	N/A	N/A	Wireless Only	N/A	Wireless Only	N/A	N/A	N/A
Picasso™	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Orange Flex Strip
Halo® iQ/ Halo 2	N/A	Cable (straight)	Cable (straight)	N/A	N/A	N/A	N/A	N/A
Muse™ iQ/ Muse	Wireless Only	N/A	Cable (straight)	Cable (straight)	Cable (straight)	Cable (straight)	Orange Flex Strip	Orange Flex Strip
SoundLens® Synergy® iQ/ SoundLens® Synergy®	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Orange Flex Strip

ADDITIONAL NOTES

- Inspire® will display any discovered non-wireless and wireless hearing aids. Non-wireless hearing aids will display first and will be automatically checked prior to initiating the fitting session. Detected wireless hearing aids must be manually selected prior to initiating the fitting session.
- VC position does not matter as the hearing aid volume control is disabled during programming.
- Starkey uses purple programming cables. Flex Strips and all boots will connect to the cables.
- Supported wired programmers: HiPro, USB HiPro, HiPro 2, NOAHLink, SpeedPort™, SpeedPort 2.
- Wireless programming requires use of either the SurfLink, NOAHLink, NOAHLink Wireless, Starkey Hearing Technologies 2.4 GHz Programmer or TruLink Programmer.

- Product can be programmed using SurfLink® Programmer

- Product can be programmed using Starkey Hearing Technologies 2.4 GHz Programmer, TruLink Programmer, or NOAHLink Wireless.

IMPORTANT

WHEN PROGRAMMING:

- It is recommended that a battery be in the hearing aid when programming with a flex strip.

PROGRAMMING TOOLS

ITEM	PART NUM.
Purple Boot	64110-000
Black Boot	64111-000
Flex Strip	65464-000

* Pictures shown are the most recent product casing.

**The RIC R is only compatible with Muse iQ.

Boots & Cables | Programming connection


Livio AI/Livio



RIC 312

BTE 13

Picasso



ITE

ITC

CIC

IIC

Muse iQ/Muse



RIC R**

RIC 312

Micro RIC 312

Standard BTE 13

Power Plus BTE 13


Mini BTE 312

ITE

ITC

CIC


Halo iQ/Halo 2



RIC 312

RIC 13

SoundLens Synergy iQ/SoundLens Synergy



IIC