

# Power BTE

Behind-The-Ear

E Series 3



#### **Sweep™ Technology**

An industry first that replaces difficult-to-adjust buttons and dials with an innovative control surface—so patients can make volume and memory adjustments with the sweep or touch of a finger

#### **Feedback Canceller**

Virtually eliminates annoying feedback

#### **Environmental Adaptation**

Continuously scans the environment and adapts appropriately for Quiet and Noise

#### **Dynamic Directionality**

Automatically adapts to ensure optimal performance in all listening situations

#### **Tonal Indicators**

Unique tones for memory, low battery, etc.


#### **4 Memories Standard**

#### **Induction Coil**


#### **Automatic Coil**

# E Series Power BTE ANSI/IEC Data


| Measurement   | Earhook                 |                    | Thin Tube<br>(Size 3+, Occluded) |                    |
|---|-------------------------|--------------------|----------------------------------|--------------------|
|   | ANSI/IEC<br>2cc Coupler | IEC OES<br>Coupler | ANSI/IEC<br>2cc Coupler          | IEC OES<br>Coupler |
| Peak OSPL90 (dB SPL)                                      | 133                     | 138                | 129                              | 132                |
| HFA OSPL90 (dB SPL)                                       | 127                     | NA                 | 113                              | NA                 |
| RTF OSPL90 (dB SPL)                                       | NA                      | 131                | NA                               | 118                |
| Peak Gain (dB)  | 70                      | 77                 | 72                               | 75                 |
| HFA Full-On Gain (dB)                                     | 66                      | NA                 | 52                               | NA                 |
| RTF Full-On Gain (dB)                                     | NA                      | 74                 | NA                               | 61                 |
| <b>Frequency Range (Hz)</b>                               | 100 - 6000              | 100 - 6000         | 100 - 5800                       | 100 - 5800         |
| Reference Test Frequency (kHz)                            | NA                      | 1.6                | NA                               | 1.6                |
| HFA Frequencies (kHz)                                     | 1.0, 1.6, 2.5           | NA                 | 1.0, 1.6, 2.5                    | NA                 |
| Reference Test Gain (dB)                                  | 50                      | 56                 | 36                               | 43                 |
| <b>Harmonic Distortion</b>                                |                         |                    |                                  |                    |
| 500 Hz (%)  | <4                      | <4                 | <1                               | <1                 |
| 800 Hz (%)  | <1                      | <1                 | <1                               | <1                 |
| 1600 Hz (%)   | <1                      | <1                 | <1                               | <1                 |
| <b>Attack and Release Time<br/>(ANSI/IEC) – Test Mode</b> |                         |                    |                                  |                    |
| Attack Time (ms)  | 22                      | 5                  | 22                               | 5                  |
| Release Time 0.1s (ms)                                    | 5-150                   | 5-250              | 5-150                            | 5-250              |
| Release Time 2.0s (ms)                                    | 5-150                   | 5-250              | 5-150                            | 5-250              |
| <b>Induction Coil Sensitivity</b>                         |                         |                    |                                  |                    |
| HFA SPLITS (ANSI) (dB SPL)                                | 112                     | NA                 | 96                               | NA                 |
| MASL (IEC) (dB SPL)                                       | NA                      | 102                | NA                               | 91                 |
| ANSI/IEC Battery Current (mA)                             | 1.8                     | 1.8                | 1.8                              | 1.8                |
| Idle Current (mA)   | 1.5                     | 1.5                | 1.5                              | 1.5                |
| <b>Estimated Battery Life<br/>for 16-Hour Day</b>         |                         |                    |                                  |                    |
| 13 Zinc Air (days)  | 9-12                    | 9-12               | 9-12                             | 9-12               |



OSPL90 (blue) and Full-On Gain (gray) curves for the E Series 3 Power BTE with Earhook.



OSPL90 (blue) and Full-On Gain (gray) curves for the E Series 3 Power BTE with Thin Tube.



E Series 3 Power BTE fitting range.

## Measurement Conditions and Recommendations

The data for E Series are obtained and performance is expressed according to ANSI S3.22 (2003), ANSI S3.7 (1995), ANSI C63.19 (2007), IEC 60118-7 (2005), IEC 60711 (1981), DIN 45605 (1989), and IEC 60118-0 (1983) with Amendment 1 (1994-01). The Starkey proprietary Real Time Analyzer and the Starkey Automated Design Verification Test System (SADVTS) comprise the basic test equipment. Data may be subject to change with product refinement.

Because of the adaptive signal processing capabilities of E Series hearing instruments, the hearing instrument must be set to test mode to compare the actual performance of the hearing instrument with these specifications. E Series hearing instruments may be set to test mode with Inspire® by reading the hearing aid and selecting the "Hearing Aid Test" screen from the menu on the left side of the Inspire window, then selecting the "Full on Gain" button.