

Starkey Hearing Technologies TV Streamer




Pairing

Starkey® Hearing Technologies TV Streamer must be paired to interact with Livio™ AI and Livio hearing aids. Pairing occurs once, then the hearing aids and TV streamer will automatically recognize each other. The TV streamer can be paired within Inspire® or outside of Inspire. Livio AI and Livio hearing aids can be paired with up to five TV streamers.


Pairing within Inspire

- 1 Launch Inspire, connect and detect hearing aids.
- 2 Select “Accessories” on the left navigation bar.
- 3 Ensure the TV streamer is powered on.
- 4 Place the hearing aids within 6 inches of the TV streamer.
- 5 Place the TV streamer in pairing mode with a short press of the pairing button. The LED will rapidly blink blue.
- 6 Select “Pair” to put the hearing aids into pairing mode (A).
- 7 Pairing is successful when the TV streamer image displays in the accessory screen (B).
- 8 Use the drop-down menu below the streaming accessories box to set “Hearing Aid Mic While Streaming” (C).
Select the information icon for pairing instructions if needed (D).



Pairing outside of Inspire

- 1 Power off any mobile phone or other accessories that were previously paired with the hearing aids.
- 2 Ensure the TV streamer is powered on.
- 3 Open and close the battery doors of the hearing aids.
- 4 Place the hearing aids within 6 inches of the TV streamer.
- 5 Press the Bluetooth® button on the TV streamer to place in pairing mode. The LED will rapidly blink blue.
- 6 Once the LED blinks green for each hearing aid, both hearing aids are paired. The TV streamer will recognize if the fitting is monaural or binaural and pair one or two hearing aids accordingly.




Basic Setup



- Ⓐ Power Supply Jack (micro USB) Ⓑ Toslink Jack Ⓒ Analog 3.5mm Jack Ⓓ Mono/Stereo Switch

- 1 Remove the plastic caps from each end of the TOSLINK cable.
- 2 Plug the TOSLINK cable into the Audio OUT jack of the electronic audio source. If there is no TOSLINK audio OUT jack, use a 3.5mm to RCA cable.
- 3 Power on the electronic audio source.
- 4 Plug the other end of the TOSLINK cable into the rear of the TV streamer Ⓑ.
- 5 Insert the power adapter microUSB connector into the rear of the TV streamer Ⓐ.
- 6 Plug the power adapter into the wall socket.
- 7 Turn on the TV streamer. A solid green LED indicates power is on.
- 8 A solid blue LED indicates the TV streamer has detected audio and is streaming.



Use

An unlimited number of Livio AI and/or Livio hearing aid users can connect to the TV streamer at one time, allowing them to share in the same listening experience.

The TV streamer supports both analog and digital input and can be successfully connected to many electronic audio sources including televisions, computers, MP3 players, DVD players, cable boxes, and Amazon® Echo and Echo Dot (3.5mm Jack).

Volume

Streamed audio volume level can be adjusted via the hearing aid user control, TV streamer, remote control or Thrive™ Hearing app.¹

Audio Streaming

Audio streaming may be started/stopped in several ways:

- Hearing Aid User Control: Push & Hold or Double Tap²
- Remote Control: Favorite Button
- Thrive Hearing app

Audio streaming will stop when:

- The TV device is powered off.
- A phone call via direct mobile phone connection or connection with Remote Microphone + is received. After the phone call concludes, all streaming stops. If desired, the TV stream will need to be manually restarted.
- If the user walks out of range of the TV stream (Up to 45 feet or 15 meters) streaming stops and will need to be manually restarted.

¹Accessories only pair to hearing aid (not to iPhone/app). ²Double Tap feature only available on Livio AI products.

Livio AI and Thrive are compatible with iPhone X, iPhone 8 Plus, iPhone 8, iPhone 7 Plus, iPhone 7, iPhone SE, iPhone 6s Plus, iPhone 6s, iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, iPad Pro, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 4, iPad mini 3, iPad mini with Retina display, iPad mini, iPod touch (5th generation) and Apple Watch. Use of Apple Watch requires iPhone 5 or later.

Apple, the Apple logo, iPhone, iPod touch, Apple Watch, App Store and Siri are trademarks of Apple, Inc., registered in the U.S. and other countries.