

QuickTIPS



Hearing Is Our Concern™

NAVIGATING THE NOISE CONTROL SCREEN

PRODUCT AVAILABILITY: Z Series, Halo, 3 Series, Xino, Xino Classic*, Wi Series,
X Series: i110/110, i90/90, i70/70

*Limited feature set in Xino Classic


SETUP

Noise control features are product, technology tier, and/or style dependent. Higher levels of technology provide more options and assistance for the patient. Default values are research-driven based on patient performance and preference.

Noise Control settings can be configured on a per-memory basis.

1. Launch **Inspire 2014.3** or higher.
2. Connect and detect hearing aids to the fitting computer via the appropriate programming accessories. Refer to the **Boots and Cables Programming Guide** for additional information regarding connecting the hearing aids.
3. Select **New Session** on the **Get Started** Screen, then select **Start**.
NOTE: Auto Path will default to launching automatically for hearing aids that have not been previous read by Inspire.
4. Select **Noise Control** on the left navigation bar.
5. Select the desired memory from the **Memory tabs** on the top right corner of the Noise Control screen.
6. If changes are desired for all memories, select the **Link Memories** icon to the left of the memory tabs to link the memories together prior to making changes.


- » Noise Control for Z Series now includes Acuity Directionality to provide consistent, high-fidelity sound that helps patients listen to speech in even the noisiest environments.
- » Z Series now offers more aggressive wind adaption than previous hearing aids.



Noise Control Screen

VOICE iQ2

- » To make Voice iQ2 more aggressive (more gain adaptation); select a radio button corresponding to more activity.
- » To make Voice iQ2 less aggressive (less gain adaptation); select a radio button corresponding to less activity.



Voice iQ2

QUIET

- » To make Quiet more aggressive (greater expansion); select a radio button corresponding to more activity.
- » To make Quiet less aggressive (less expansion); select a radio button corresponding to less activity.



Quiet

Fast-acting noise management and speech preservation system intended to help patients hear effectively in noise and reduce listening effort.

Expansion algorithm designed to provide comfort for low level noise.

WIND

- » To make Wind more aggressive (more gain adaptation); select a radio button corresponding to more activity.
- » To make Wind less aggressive (less gain adaptation); select a radio button corresponding to less activity.
- » To alter Wind Time Constant (speed of engagement):
 - Select **Wind** feature name
 - Select **Slower** or **Faster** as desired

NOTE: Alteration of Wind Time Constant is available in all i110/110 and i90/90 levels of technology (with the exception of Xino Classic).

Noise management system designed to provide comfort for wind noise when turbulence is detected over the hearing aid microphones.




MACHINE NOISE

- » To make Machine Noise more aggressive (more gain adaptation); select a radio button corresponding to more activity.
- » To make Machine Noise less aggressive (less gain adaptation); select a radio button corresponding to less activity.
- » To alter Machine Noise Time Constant (speed of engagement):
 - Select **Machine Noise** feature name
 - Select **Slower** or **Faster** as desired



NOTE: Alteration of Machine Noise Time Constant is available in all i110/110 and i90/90 levels of technology (with the exception of Xino Classic).

Noise management system designed to provide comfort for loud, steady state noise.



DIRECTIONALITY (MICROPHONE MODE)

1. Select desired **Directionality** setting from drop-down menu:
 - » **Omni** (Omnidirectional): Fixed Response; amplifies sound from all directions equally.
 - » **Directional** (Fixed directional): Fixed Response; amplifies sound from in front of the listener more than from behind via a hypercardioid polar plot.
 - » **Dynamic** (Omni and Directional): Automatic switching; engages in omni or fixed directional based on environment.
 - » **Adaptive** (Adaptive Directional): Automatic switching; Acuity Directionality with null steering and Speech ID to protect speech at all angles around the listener.
2. Select **Directionality** feature name for Advanced Directionality settings:
 - » Select **Less Often** or **More Often** under Switching as desired.
 - » Select **Slower** or **Faster** under Directionality Time Constants as desired.



Advanced Directionality Settings

DIRECTIONALITY PLUS

1. Select **Directionality Plus** to enable.
 - » Defaults OFF.
 - » Available in Z Series, 3 Series wireless, and Xino wireless: i110, i90, and i70 products.

- » A more aggressive noise management strategy, providing additional gain adaptation for noise when speech is present.
- » Directional response must be engaged via Dynamic or Adaptive Directionality.

BINAURAL SPATIAL MAPPING

1. Deselect **Binaural Spatial Mapping 3D** to disable for Wind, Machine Noise, and/or Directionality.
 - » Defaults ON.
 - » Available in all i110, i90, and i70 wireless products (with the exceptions of Wi Series i70 and the Halo product family).

- » Ear-to-ear collaboration between two wireless hearing aids in a binaural fit to improve comfort and provide the best SNR possible for the listener.
- » Engages Wind, Machine Noise, and Directionality.
- » Tier-dependent

GO LIVE

1. Select **Go Live** to view Noise Control settings in real time.
2. Simulate the desired environment by selecting audio files from the **Media Player** or other source.
3. Make necessary adjustments.
4. Deselect **Go Live** to disengage the tool.

A tool for assisting the professional in adjusting the noise control settings by providing a real-time display of noise control engagement.



Go Live