
TRULINK™
PROGRAMMER

INSTALLATION INSTRUCTIONS

TABLE OF CONTENTS

Description	4
Installation	4
Operation	5
Safety Precautions	6
Regulatory & Service Information	6
Important Safety and Usage Information	7
Regulatory Notices	9

This document contains instructions for connecting TruLink Programmer into your laptop or desktop PC via a USB connection. TruLink Programmer requires that Inspire® 2014.2 fitting software or higher is installed.

DESCRIPTION

TruLink Programmer is a USB hearing aid programmer functioning as the interface between wireless hearing aids and computer fitting software. It consists of a USB docking station and TruLink programmer dongle.

INSTALLATION

NOTE: Screen may vary depending on the Windows® version.

1. Confirm that Inspire 2014.2 or higher is installed. (TruLink Programmer is not supported on earlier versions of Inspire fitting software.)
2. Insert the USB cable of the docking station into an open USB port on your PC.

NOTE: It is strongly recommended to use the programmer as it is assembled (dongle plus USB docking station). Plugging the dongle directly into the computer is not recommended.

3. TruLink Programmer is designed to be “plug and play” and does not require a CD or a restart of the PC. Device drivers will install automatically.

OPERATION

Once installed, TruLink Programmer is ready to use.

To use:

1. Verify TruLink Programmer is connected to the PC
2. Launch Inspire 2014.2 or higher
3. Search for TruLink Programmer, if needed, by following the steps below:
 - *Select Tools*
 - *Select Programmers*
 - *Click Search*
 - *Select TruLink Programmer to program the hearing aids*

SAFETY PRECAUTIONS

Please read the following safety precautions before using TruLink Programmer for the first time.

Markings and Symbols

The following markings and symbols are used in these safety precautions and/or on the device label.



Within Europe, dispose of TruLink Programmer according to local EU WEEE regulations.

REGULATORY AND SERVICE INFORMATION

Classification

- Classification to EN 60950-1 Class II Continuous Operation.
- Any associated computer used with TruLink Programmer must be an EN 60950 (ISO 950) approved computer.

Power Supply

- Operating Voltage: 5.0 vDC supplied from the USB port of an EN 60950 (ISO 950) approved computer.

IMPORTANT SAFETY AND USAGE INFORMATION

Storage and Shipping Conditions:

- Temperature: -40 to 85° C, -40 to 185° F
- Humidity: 30 – 90%
- No condensation

Operating Conditions:

- Temperature: -40 to 85° C, -40 to 185° F
- Humidity: 30 – 90%
- No condensation

Avoid using TruLink Programmer in areas subject to the following conditions:

- Extreme heat or cold
- Static electricity
- Large amounts of dust
- Rough handling

Rough handling or dropping TruLink Programmer can cause damage. If you think you may have damaged TruLink Programmer, or if TruLink Programmer is suffering an obvious malfunction, immediately unplug it from your computer and contact your authorized dealer to request servicing.

DO NOT continue to use a damaged TruLink Programmer.

DO NOT place heavy objects on top of TruLink Programmer.

Never insert foreign objects into the connections.

Avoid spilling liquid on TruLink Programmer.

The exterior of TruLink Programmer may be wiped with a damp cloth or sterile hand wipe. Never use paint thinner, benzene, cosmetics or other volatile agents to clean the exterior of TruLink Programmer.

Never try to take apart or modify TruLink Programmer. Doing so creates the danger of fire and electrical shock, as well as TruLink Programmer malfunction. There are no user-serviceable parts inside TruLink Programmer. All services must be performed by the authorized dealer.

REGULATORY NOTICES

FCC ID: QOQBLED112

IC: 5123A-BGBLED112

KCC: CRM-BGT-BLED112



R

003 WW A111471

FCC NOTICE

This device complies with part 15 of the FCC rules and with RSS-210. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation of the device.

NOTE: The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications to this equipment. Such modifications could void the user's authority to operate the equipment.

Hereby, Paradigm declares that TruLink Programmer is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. A copy of the Declaration of Conformity can be obtained from the below addresses.

This Class B digital apparatus complies with Canadian ICES-003.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Paradigm, Inc.
6700 Washington Ave. South
Eden Prairie, MN 55344 USA



Paradigm, Inc. European Headquarters
Wm. F. Austin House, Bramhall Technology Park
Pepper Road, Hazel Grove, Stockport SK7 5BX
United Kingdom



Waste from electronic equipment must be
handled according to local regulations.

