

# 2019 Product portfolio reference guide



## Thrive™ 2.4 GHz + NFMI<sup>1</sup> TECHNOLOGY


## STARKEY® HEARING TECHNOLOGIES ACCESSORIES

livio<sup>AI</sup>R

livio<sup>AI</sup>

livio<sup>R</sup>

livio<sup>®</sup>



The easy-to-use Thrive Hearing Control app features universal Bluetooth® connectivity to our new line of patient-focused accessories. The Thrive app is available in the app stores for both iOS and Android™ users.



Stream audio from a TV or other electronic audio source directly to Livio AI or Livio hearing aids. It offers excellent sound quality, and is easy to use.



Easily stream from a variety of audio sources. This simple-to-use remote microphone offers the most capabilities available in any accessory.



Enjoy one-on-one conversations in noisy environments with the small, easy-to-use Mini Remote Microphone. Clip it onto the clothing of the person you're talking to, or use it as a TV streamer by placing it near the speaker.



Our remote includes updated features so users can control memory and volume, mute their hearing aids or turn other features on and off.

Feature	AI Premium 2400 Healthable Technology	Premium 2400	Advanced 2000	Select 1600	Low 1200	Basic 1000
Platform: Thrive	●	●	●	●	●	●
App: Thrive Hearing Control	●	●	●	●	●	●
Healthable™ Technology: Brain and Body Tracking, Fall Alert, Heart Rate <sup>2</sup>	●					
Personal Assistance: Thrive™ Virtual Assistant	●	●				
Translate	●					
Transcribe	●					
Self Check	●	●	●	●	●	●
Tap Control	●					
Telehealth Service: Hearing Care Anywhere™	●	●	●	●	●	●

Feature	AI Premium 2400 Healthable Technology	Premium 2400	Advanced 2000	Select 1600	Low 1200	Basic 1000
Sound Imaging: Channels   Bands	24	24	20	16	12	10
Bandwidth	10 KHz	10 KHz	10 KHz	10 KHz	8 KHz	8 KHz
Speech Optimization	●	●	●	●	●	●
Music Optimization	●	●	●	●	●	●
Music Adaptation	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
E2E Music Adaptation	●	●	●	●	●	●
Ear-to-Ear Technology	●	●	●	●	●	●
E2E Wind Noise Management	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
E2E Machine Noise Adaptation	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
E2E Directionality	●	●	●	●	●	●
E2E Phone Streaming <sup>2</sup>	●	●	●	●	●	●
Sound Manager	●	●	●	●	●	●
Music	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Speech	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Speech and Noise	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Machine Noise	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Wind	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Noise	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Quiet	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4

Feature	AI Premium 2400 Healthable Technology	Premium 2400	Advanced 2000	Select 1600	Low 1200	Basic 1000
Directional Processing	●	●	●	●	●	●
Channels	24	24	20	16	12	10
Immersion	●	●	●	●	●	●
Adaptive	●	●	●	●	●	●
Dynamic	●	●	●	●	●	●
Directional	●	●	●	●	●	●
Noise Control	●	●	●	●	●	●
Noise Reduction Strength	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Spatial Speech Enhancement	●	●	●	●	●	●
Transient Noise Reduction	●	●	●	●	●	●
Feedback Management	●	●	●	●	●	●
Frequency Lowering	●	●	●	●	●	●
CROS System	●	●	●	●	●	●
Tinnitus Technology: Multiflex Tinnitus Technology	●	●	●	●	●	●
Compatibility: Starkey Hearing Technologies Accessories	●	●	●	●	●	●

<sup>1</sup>Excludes micro RIC 312

<sup>2</sup>Coming soon for Livio AI RIC R hearing aids with Custom Healthable Receiver

# Synergy® 900SYNC TECHNOLOGY

## Muse iQ<sup>R</sup>



RIC R<sup>3</sup>



Feature	Premium i2400	Advanced i2000	Select i1600
Platform: <b>Synergy</b>	●	●	●
Sound Imaging: <b>Channels   Bands</b>	24	20	16
Bandwidth: <b>Extended 10 KHz</b>	●	●	●
Speech Compression: <b>Acuity Speech Optimization</b>	●	▶	▶
Spatial Awareness: <b>Acuity Immersion<sup>4</sup></b>	●	●	●
Music Optimization: <b>Music Enhancement</b>	●	●	●
Music Adaptation <b>E2E Music Adaptation</b>	1 2 3 4	1 2 3 4	1 2 3 4

## Muse iQ



RIC 312t



micro RIC 312t



mini BTE 312



Power Plus BTE 13



Feature	Premium i2400	Advanced i2000	Select i1600
Ear-to-Ear Technology: <b>Acuity Binaural Imaging</b>	●	▶	▶
E2E Wind Noise Management	1 2 3 4	1 2 3 4	1 2 3 4
E2E Machine Noise Adaptation	1 2 3 4	1 2 3 4	1 2 3 4
E2E Directionality	●	●	●
E2E Phone Streaming	●	●	●
Environment Manager: <b>Acuity Lifescape Analyzer</b>	●	●	●
Music	1 2 3 4	1 2 3 4	1 2 3 4
Speech	1 2 3 4	1 2 3 4	1 2 3 4
Speech and Noise	1 2 3 4	1 2 3 4	1 2 3 4
Machine Noise	1 2 3 4	1 2 3 4	1 2 3 4
Wind	1 2 3 4	1 2 3 4	1 2 3 4
Noise	1 2 3 4	1 2 3 4	1 2 3 4
Acuity Quiet	1 2 3 4	1 2 3 4	1 2 3 4

## Picasso™



ITE



ITC



CIC



IIC<sup>5</sup>



Feature	Premium i2400	Advanced i2000	Select i1600
Directional Processing: <b>Acuity Immersion</b>	●	▶	▶
Channels	24	20	16
Immersion	●	●	●
Adaptive	●	▶	▶
Dynamic	●	▶	▶
Directional	●	▶	▶
Noise Control: <b>Acuity Voice</b>	●	▶	▶
Noise Reduction Strength	1 2 3 4	1 2 3 4	1 2 3 4
Feedback Management: <b>WhistleFree Feedback Cancellation</b>	●	●	●
Frequency Lowering: <b>Speech Shift</b>	●	●	●
Tinnitus Technology: <b>Multiflex Tinnitus Technology</b>	●	●	●
Compatibility: <b>SurfLink Wireless Accessories</b>	●	●	●

# SURFLINK® ACCESSORIES



### SurfLink Mini Mobile

Patients can stream calls via Bluetooth to their Muse iQ, Muse iQ R and wireless Picasso hearing aids.



### SurfLink Remote Microphone 2

Stream audio directly to Muse iQ, Muse iQ R and wireless Picasso hearing aids or use to help patients in one-on-one conversations.



### SurfLink Media 2

Just plug a TV, radio or MP3 player into SurfLink Media and rich, stereo sound is streamed directly to any SurfLink compatible hearing aids.



### SurfLink Mobile 2

Provides patients with the convenience of hands-free cell phone conversations by turning our wireless hearing aids into the phone microphone and receiver.



### SurfLink Remote

Patients can control memory and volume adjustments and even adjust their hearing aids, or easily go in and out of streaming mode.

<sup>3</sup> Also available in i1200 and i1000 technology tiers

<sup>4</sup> CIC style only

<sup>5</sup> Only available in the premium technology tier



Tinnitus Technology



Absolute Power/  
Power Products



CROS System



Rechargeable Option

For Livio AI, Livio and Thrive Hearing Control app and iOS/Android compatibility information, please visit [starkey.com/thrive-hearing](http://starkey.com/thrive-hearing).

Use of the Made for Apple badge means that an accessory has been designed to connect specifically to the Apple product(s) identified in the badge, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Apple, the Apple logo, iPhone, iPad, iPod touch, App Store and Siri are trademarks of Apple Inc., registered in the U.S. and other countries.

Google Play and Android are trademarks of Google Inc.

All trade names and trademarks are properties of their respective owners.