

# AURACARE HEARING AID OPTIMIZATION SYSTEM

Of all that you do, perhaps nothing is more important to your long-term viability than the quality of service you provide to each client. So how can you raise your level of care and distinguish yourself in today's increasingly competitive hearing care environment?



With the **AuraCare Hearing Aid Optimization System**, you have the procedures and technology necessary to perform comprehensive examinations of a hearing aid's physical and electrical condition. All at your fingertips in one convenient system.



## INTEGRITY METER

A special probe measures battery drain, transducer functionality, and the consistency of volume control taper.



## APR CHAMBER

Through the use of acute atmospheric pressure reduction, the chamber helps accelerate the removal of internal moisture.

## COMPARITOR

Lets your clients listen to a unique side-by-side sound comparison of their aid's performance with the new, Sequel™ High Performance hearing aid.



## AURAVAC

Allows you to focus a powerful pinpoint vacuum to extract obstructions and debris in microphone, receiver, and vent ports.



HEARING CARE EXCELLENCE



AURACARE HEARING AID OPTIMIZATION SYSTEM