

Harp Plus
now available
with optional
QuickSIN and
TEN Test

DIAGNOSTIC AUDIOMETER



Harp



Audiology
line

Harp

Harp is an advanced diagnostic audiometer, capable of performing fast and accurate air, bone and speech audiometric exams, including free field.

Impressive features combined with ease of use make Harp the ideal choice for a wide range of users, from private doctors to hospitals and clinics.

COMPUTER INTEGRATION

At Inventis we see computer integration as a key feature in all our devices. That is why we have developed Daisy: modular software for storing and managing all kinds of ENT and audiology examinations.

Simply connect Harp to a USB port on your computer (no SW driver needed), and you can transfer exam results into patient folders in the Daisy database, complete with, e.g., video otoscope or vestibular session files. All diagnostic data is therefore accessible in just a few clicks.

A NOAH module for Harp is also available.

TECHNICAL DATA

- Intensity AC: -10 to 120 dB HL. BC: -10 to 80 dB HL
- Frequency range... 125-8000 Hz
- Input Tone, warble tone, external source (tape, CD) 1-2, mic., flash memory (only in Plus model)
- Masking..... NBN, WN, SN
- Output..... AC L+R, BC, FF L+R, insert masking
- Talk over..... Built-in or external microphone
- Monitor..... Built-in speaker or ext. headset or speaker
- Display..... Colour graphical display, 95 x 54 mm, 480 x 272 pixels
- Test battery Speech audiometry, Autothreshold, ABLB.
Only in Plus model: SISI, DLI, Tone decay, Stenger, Master Hearing Aid
- Data storage..... Internal memory for 100 patients
- Printer..... Optional integrated thermal printer
- PC interface..... USB
- Compatible SW..... Inventis Daisy or Noah with dedicated module
- Dimensions 32 x 32 x 9 cm / 12.6 x 12.6 x 3.5 in. (LxWxH)
- Weight 1.8 Kg / 4.0 lbs
- Standards..... Audiometer: EN 60645 -1 / ANSI S3.6, Type 2
Speech: EN60645-2 / ANSI S3.6 type A-E
Safety: EN 60601-1, Class 2, Type BF
EMC: EN 60601-1-2
- Classification..... Class IIa (MDD 93/42)

THE MODELS

Harp is available in two configurations to suit different user needs: Basic and Plus.

THE TESTS

Besides standard audiometry exams, Harp also performs Hughson-Westlake automated audiometry, speech audiometry and ABLB exams. The "Plus" model also performs SISI, Stenger (pure tone and speech), DLI and Tone Decay tests and includes a 2 channels Master Hearing Aid.

EMBEDDED SPEECH MATERIAL

No more CD player headaches! Harp Plus features an integrated 2 GB flash memory for storing your speech material. Dedicated buttons allow you to access word lists with ease. The screen also displays the word currently being pronounced.

PORTABLE SYSTEM

Despite its top-level features, Harp is ultra-compact and easily transportable. It even comes with a dedicated carrying case (optional).



THE DISPLAY

Harp features a large, bright graphic colour display, on which you can view the complete audiogram during the examination.



6700 Washington Ave S
Eden Prairie, MN 55344

800.733.2636
800.328.8602
952.941.6401

© Starkey Hearing Technologies.
All Rights Reserved.
BROC0319-00-EE-XX 6/12

Harp is developed by:

INVENTIS s.r.l.
CORSO STATI UNITI, 1/3
35127 PADOVA - ITALY
PHONE: +39.049.8962 844
FAX: +39.049.8966 343
info@inventis.it
www.inventis.it

Follow us on Twitter
<http://twitter.com/inventissrl>

 The Inventis Quality System complies with ISO 9001 and ISO 13485 standards.

Inventis® is a registered trademark of Inventis s.r.l.