



STARMED PROFESSIONAL FITTING STATION

Your Commitment To Optimal Hearing



TECHNICAL SPECIFICATIONS



STARMED Professional Fitting Station

The specially designed Professional High Speed Oto-Lathe is built for ease of use and rugged dependability. It can help you quickly and precisely modify any style of hearing aid or earmold for improved fit and comfort. A special industrial grade diamond bit provides smooth, injury free plastic reduction. The QuickChange Chuck makes for fast and simple bit replacement. Additional burrs for enlarging or modifying standard vents plus special flexible bits for CIC vent modification are also included. The Adjustable Work Light and Magnification Lens can be positioned up or down to ease the strain of any precision modification. Each Buffing Wheel is pre-treated and compounded for easy scratch removal and high polish, helping restore your patients' hearing aids to like-new condition. The wheel is surrounded by the Padded Safety Bin to prevent damage to the hearing aid. The StarMed EarLab makes UV shell repairs and modifications quick and easy. The proprietary Ultra Shell Material hardens instantly in the UV Curing Chamber to produce a professional, long lasting repair. The Ultra Shell Process may also be used to add plastic to improve comfort or reduce feedback. The multi-purpose StarMist material can handle a wide variety of repair and modifications, including the patching and repair of dental acrylic shells, as well as forming and attaching CIC removal handles. The Professional Fitting Station also contains an integrated Pull-Out Work Shelf to accommodate more time-consuming repairs or modifications. Also included in the drawers are all of the tools and modification kits necessary for the professional fitting and service of hearing instruments. Plus, the station's large, heavy duty casters make for easy mobility on all types of floor surfaces.