

INTRODUCING

Halo™ BTE 13

Our Newest **Made for iPhone®**
Hearing Aid.



The most personalized
hearing experience
continues.



Hearing Is Our Concern™

Connect to
HEARING.

Connect to
HEALTH.

Connect to
LIFE.



Connect to Halo.

Connect to LIFE.



*Help your
patients
hear better,
live better.*


We are proud to introduce the new **Halo BTE 13** and the latest advances to our Made for iPhone® Hearing Aids. Now you'll be able to fit even more patients with the most revolutionary and personalized form of wireless connectivity available today.

With the addition of the BTE 13 –the Halo family continues to break new ground with overall performance improvements and sound enhancements, including our new evidence-based, patient-focused breakthrough operating system – **BluWave® 4.0**.

*Halo BTE 13
features up to
70dB gain for
patients with
more severe
hearing loss.*

**Keep your patients
active and connected
with BluWave 4.0.**

At the heart of all Halo products is our benchmark operating system – BluWave 4.0. This new operating system's signal processing power enables us to continually optimize audibility, comfort and sound quality. The result is **Starkey Superior Sound™**.



Inspire® 2015.1 is required
to fit the new Halo BTE
and BluWave 4.0 updates.



*Your patients' overall health can be
directly linked to their hearing health*

Connect to the TruLink Hearing Control App

We've also updated our patient-preferred **TruLink™ Hearing Control app** to include features that enhance streaming and hearing in noise. With **Stream Boost**, patients can have an automatic setting that boosts the incoming media stream for enhanced audio performance. **Comfort Boost** allows patients to aggressively reduce noise while optimizing sound quality in noisy environments.

We have introduced all these revolutionary technology updates so Halo will be a hearing aid that helps your patients **hear better, live better** and have a healthier life.

Push Button – BTE features new push button feature for on-board user control.

Halo and TruLink are compatible with iPhone 6, iPhone 6 Plus, iPhone 5s, iPhone 5c, iPhone 5, iPhone 4s, iPad Air, iPad Air 2, iPad (4th generation), iPad mini with Retina display, iPad mini, iPad mini 3 and iPod touch (5th generation).



Patients want ... Halo delivers

A More Natural Listening Experience

Starkey's **Binaural Spatial Mapping** delivers maximum listening comfort for louder sounds while providing a more natural, overall sound. **Voice iQ™2** is a noise management and speech preservation system intended to let patients hear effectively in noise.

Superior Audibility in Noisy Environments

Thanks to **Acuity Directionality™**, our new multichannel adaptive directionality solution, speech protection is improved in difficult listening situations. Starkey's proprietary **Speech ID** detects speech even when it occurs behind the patient. The combination of these two innovative technologies results in consistent, high-fidelity sound that helps patients listen to speech and participate in conversations in even the noisiest environments. The system's Micro-Electronic Mechanical Systems (MEMs) microphones are the same as found in today's iPhone.

Distortion-Free Comfort in Loud Situations

Our **ISO-Clear Compression** delivers a more personalized listening experience by making soft speech sounds audible without affecting louder inputs. Loud sounds are discreetly managed while maintaining the quality of lower-level sounds.

Say Farewell to Feedback

Starkey's best-in-class feedback cancellation system, **PureWave Feedback Eliminator** has the ability to manipulate how aggressively it performs, so feedback elimination can be personalized to the specific product and patient. It's designed to identify, classify and eliminate annoying whistling and buzzing in all environments, so your patients experience feedback-free listening.

Laughs at Water and Wax

Halo features Starkey's proprietary moisture protection, **HydraShield®2**. When combined with Halos' unibody construction, this nano-coating delivers excellent protection from water, wax and oils.

Enjoy Every Conversation

Engineered to enhance real-time audibility by identifying high-frequency speech cues and replicating them in lower frequencies, **Spectral iQ** maintains frequency relationships retaining sound quality while providing clarity for high-frequency speech.

Join the Crowd With Active Noise Reduction

Comfort Boost allows patients to aggressively reduce noise on-demand while optimizing sound quality in environments like noisy restaurants or crowded rooms.

Make Control Adjustments a Thing of the Past

Auto Experience Manager provides gradual increased gain settings over a period of time to help first-time hearing aid patients transition to Halo.



Hearing Is Our Concern™

**6700 WASHINGTON AVENUE SOUTH
EDEN PRAIRIE, MN 55344
1.800.328.8602**

StarkeyPro.com

"Made for iPod," "Made for iPhone," and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, or iPad may affect wireless performance.

Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.