

# IGNITE™

.....  
**WIRELESS HEARING  
PRODUCTS**



Hearing Is Our Concern™

# YOUR JOURNEY TO BETTER HEARING STARTS HERE

**ignite**<sup>™</sup>  
WIRELESS

When it comes to wanting small, comfortable and great sounding hearing aids, more people are turning to Starkey, the only American-owned and operated manufacturer of advanced hearing solutions.

And for many first-time wearers, Ignite Wireless hearing aids are the perfect choice.



Ignite Wireless is Starkey's full line of value-priced hearing aids that sound great and perform dependably in a variety of normal listening situations. Engineered by Starkey, America's premier hearing aid provider, all Ignite Wireless hearing aids are designed to deliver superior:

**Feedback Cancellation** – greatly reduces annoying buzzing and whistling caused by things like talking on the phone or being hugged.

**Noise Management** – automatically recognizes and reduces noise for more comfortable listening in loud environments.

**Directional Settings** – designed to help you understand the speech in front of you while filtering out unwanted background noise.

# HEAR ALL THAT LIFE HAS TO OFFER

---

Combined with our innovative SurfLink® accessories, Ignite Wireless hearing aids can also make it easier to enjoy TV, music, cell phone conversations and more — no wires attached!

## ↘ SURFLINK MEDIA

This set-and-forget media streaming solution connects to TVs, MP3 players and more, to wirelessly stream audio directly to your Ignite Wireless hearing aids. This allows you to hear the volume at the level you want, while everyone else in the room can listen at the volume they want.



## ↘ SURFLINK REMOTE

This handy remote lets you adjust and control your Ignite Wireless hearing aids with the push of a button.





## ↘ SURFLINK MOBILE

Enjoy hands-free phone conversations with this one-of-a-kind device. That's because SurfLink Mobile turns your Ignite Wireless hearing aids into both your cell phone microphone and receiver:



**Microphone** — Your hearing aids can pick up your voice and send it to the person on the other end of the phone.

**Receiver** — SurfLink Mobile streams the voice on the other end of the phone directly to both of your hearing aids

But that's just the beginning. SurfLink Mobile is not only a cell phone transmitter, it's also an assistive listening device, media streamer and hearing aid remote all rolled into one.

# WHICH IGNITE WIRELESS IS RIGHT FOR YOU?

Ignite Wireless hearing aids are available in a full range of styles, from In-The-Ear (ITE), In-The-Canal (ITC), Half Shell (HS), and Starkey's renowned Completely-In-Canal (CIC) to standard Receiver-In-Canal (RIC) devices.





ignite™  
WIRELESS





Hearing Is Our Concern™

[www.Starkey.com](http://www.Starkey.com)

---

6700 WASHINGTON AVE. S.  
EDEN PRAIRIE, MN 55344  
1.800.328.8602



Individual results may vary.