

TREAT TINNITUS

THE PERSONALIZED WAY



XINO™
TINNITUS

INTRODUCING A MORE FLEXIBLE WAY TO TREAT TINNITUS

xino[™]
TINNITUS



Few ailments are as subjective as tinnitus. From what it sounds like to when it's most annoying to how debilitating it can be, tinnitus affects every person differently.

New Xino Tinnitus tackles the subjective nature of tinnitus head-on — all in a nearly invisible RIC 10 device.

Using our patent-pending Multiflex Tinnitus Technology, Xino Tinnitus creates a patient-preferred sound stimulus designed to deliver the one thing every sufferer wants — relief. It makes tinnitus manageable for patients who have been frustrated by their inability to find consistent, personalized relief — and gives hearing professionals an easy-to-use tinnitus management option.

PATIENTS WANT:



XINO TINNITUS DELIVERS:

**Relief for ringing
in the ears**

A discreet RIC device designed to help patients with tinnitus

**A way to customize
relief to match their
unique tinnitus**

Multiflex Tinnitus Technology, which enables superior flexibility in fine-tuning a relieving noise stimulus

**A tinnitus solution
and hearing aid
rolled into one**

All the performance features of our high-end hearing aids as well as Multiflex Tinnitus Technology

**A fully configurable
solution**

A device that can be configured (on a per-memory basis) as hearing aid only, hearing aid + tinnitus stimulus, or tinnitus stimulus only

**Quick and easy,
device-level control**

On-the-fly volume and memory adjustments to either tinnitus stimulus or the hearing aid

**Compatibility with
multiple tinnitus
therapy protocols**

Seamless integration with any sound therapy, including retraining or masking

**A way to deal with
and manage the
subjective nature
of tinnitus**

SoundPoint for Tinnitus, which helps the patient identify and communicate sound preferences at the fitting stage

• Xino Tinnitus combines advanced
• hearing aid technology with
• personalized tinnitus relief.
•

Because tinnitus is different for every patient, we designed Xino's Multiflex Tinnitus Technology with personalization and flexibility in mind. We also designed it to seamlessly integrate into most tinnitus treatment strategies.

Combined with Inspire® 2013.3, Xino Tinnitus gives you unprecedented ability to uniquely generate a sound stimulus that provides the greatest relief:

- First, our "best-fit" tinnitus algorithm intelligently sets the tinnitus stimulus according to the patient's audiogram.
- Next, 16 bands of frequency adjustments and our flexible modulation algorithm make it easy to shape the stimulus to fit your patient's unique tinnitus experience.
- Finally, SoundPoint Tinnitus allows patients to become an active partner in fine-tuning the stimulus to their precise preference, which helps ensure greater success.



Equipped with industry leading technology and our new, highly customizable Multiflex Tinnitus Technology, Xino Tinnitus offers a solution for the millions of patients who have ringing in their ears.



Customize relief to your patient's unique tinnitus experience

Key Xino Tinnitus features:

- Tinnitus therapy + hearing aid = two devices in one
 - Configured memories (tinnitus, tinnitus + hearing aid, hearing aid)
- Nearly invisible device – RIC 10
- Starts with a best fit and then allows for flexible adjustments in 16 bands
 - Amplification operates independently of tinnitus stimulus
- Ability to adjust modulation rate of stimulus
- Flexible device-level control of volume and memory
- Can adhere to any fitting philosophy, including retraining therapy and masking protocols
- SoundPoint Tinnitus fitting tool allows patients to identify the most preferable stimulus, which helps ensure satisfaction

XINO TINNITUS
RIC 10

SoundPoint Tinnitus
increases your patient's
chance at success



WE'RE MAKING
BETTER
HEARING
MORE PERSONAL
THAN EVER

www.StarkeyPro.com



Hearing Is Our Concern™