


**FIT A LOT
IN SOMETHING
LITTLE**



**XINO™
WIRELESS**

STILL TINY.
STILL MIGHTY.
NOW WIRELESS!



 Give more patients the wireless functionality they want in the small form they demand with new Xino Wireless.

The newest member of our popular wireless lineup, Xino Wireless is a newly designed micro RIC 312 that makes it easy for even more of your patients to enjoy the features and freedom that our patented IRIS™ Technology enables.

PATIENTS WANT:



XINO WIRELESS DELIVERS:

Small hearing aids

An incredibly small RIC 312

A way to hear better in noise

Two technologies—Voice iQ² and InVision Directionality—which, together, tackle the issue of hearing speech in noise

Hands-free cell phone conversations

SurfLink[®] Mobile, the only cell phone solution that uses JustTalk to turn the patient's wireless hearing aids into both phone receiver *and* microphone

Wireless media connectivity — without any body-worn relay devices

Dongle-free wireless streaming of TV, computers, MP3 players and more, thanks to SurfLink Media and SurfLink Mobile

An easy and comfortable fitting process

The ultimate speed and flexibility in direct-to-hearing aid wireless programming via SurfLink Programmer

More natural 3D-like sound

Binaural Spatial Mapping, Starkey's proprietary ear-to-ear protocol that optimizes SNR settings based on inputs to both ears

Improved high frequency audibility

Spectral iQ, our patented technology that intelligently identifies high-frequency speech cues and replicates them in lower frequencies

A more personalized listening experience

ISO Compression, which allows for independent control of soft, moderate and loud level speech


No buzzing or whistling

The industry benchmark for feedback elimination—PureWave Feedback Eliminator— with up to 28dB of added stable gain

Reliable, no-hassle hearing aids

HydraShield^{®2}, our water and wax-repellent nano-coating; a new, ergonomic switch; and a redesigned mic cover that resists clogging



 IRIS Technology, our proprietary wireless platform, is preferred by patients and professionals thanks to its simplicity, clarity and compatibility with our innovative lineup of SurfLink accessories.

New Xino Wireless delivers everything people love about IRIS Technology — including Binaural Spatial Mapping, no body-worn relay device, hands-free wireless phone streaming, media streaming and programming, synchronized user controls and more — in a style that's comfortable, practical and desirable: a micro RIC 312.

.....

**Now patients can get
it all *and* get it small**

.....



Key Xino Wireless features:

- An incredibly small, wireless micro RIC 312
- Flexible and powerful — can fit up to 70 gain
- Compatible with all SurfLink accessories, including SurfLink Mobile
- Fully configurable, ergonomic switch enables easy volume and memory adjustments
- Loaded with Starkey's industry-leading performance features, including:
 - **Binaural Spatial Mapping** — our ear-to-ear protocol that delivers a more natural 3-D like listening experience
 - **Voice iQ²** — our superior performing noise reduction and speech preservation system
 - **InVision Directionality** — directional system that is designed to enhance speech audibility, even in complex listening environments
 - **Spectral iQ** — the industry's most intelligent way to enhance high-frequency audibility
 - **PureWave Feedback Eliminator** — adds up to 28dB of stable gain



..... Xino Wireless
is our smallest
RIC 312 ever

WE'RE MAKING
BETTER
HEARING
MORE PERSONAL
THAN EVER

www.StarkeyPro.com



Hearing Is Our Concern™